

Product Name :
Scrap Bailing Press

Product Code :
NLAB-MACHINELABWK14006



Description :

Scrap Bailing Press

Technical Specification :

Scrap Bailing Press .

RIPPLE COMPRESSION

Baling Chamber (WxHxL)

Bale Size 8-60 (WxHxL)

Electric Motor 1440 RPM

PLUNGER PUMP

Front Lid

Front Main Cylinder

Front Side Cylinder

Front Ejector Cylinder

Working Pressure

Weight of (Unit) (in MS)

000 Selection (Approx)
00000 (Approx) Machine
000 Tank Capacity NCLO-68

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003