

**Product Name :**  
Sedimentometry

**Product Code :**  
NLAB-EDUCATIONTVET13017

**Description :**

Sedimentometry

**Technical Specification :**

particle size - densimeter method according to NF standard

The set includes:

Hydrometer 0.995 to 1.030 g / ml

Standardized cylinder Ø85 mm +/- 5 mm of 2500 ml, graduated at 500, 1500 and 2000 ml

Thermometer 0 to + 50 ° C, div. 0.5 ° C

Glass tank 600x300x380 mm

1000W thermostat

Manual agitator, l = 600 mm

Mechanical agitator, 10,000 rpm with stainless steel container

250 ml beaker

Sodium hexametaphosphate, 1000g

Distilled water

Amber glass reagent bottle, 1000 ml

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

## **Naugralabequipments**

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)