

Product Name :
Seed Dryer Cabinet

Product Code :
NLABMEDICAL-M3170003

Description :

Seed Dryer Cabinet India, Seed Dryer Cabinet manufacturer India, suppliers, Medical Seed Dryer Cabinet India, Hospital Seed Dryer Cabinet India

Technical Specification :

Double walled construction with gasket double wall door.

Inner chamber made of anodized aluminum or stainless steel sheet and outer made of mild steel sheet duly finished in powder coating.

Two adjustable 90 x 89 x 5 cms.

Trays are provided.

A blower fan has been provided to maintain inner temperature uniformly.

Temperature range 5°C above ambient to 60° C ± 2° C controlled by a digital temp.

Indicator & controller.

The cabinet is mounted on castor wheels.

Optional:

Digital temp. Indicator-cum-controller with computer data logger

Chamber size:

48" x 36"

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com