

Product Name :
Seed Germinator Double Chamber

Product Code :
NLABMEDICAL-M3190003

Description :

Seed Germinator Double Chamber India, Seed Germinator Double Chamber manufacturer India, suppliers, Medical Seed Germinator Double Chamber India, Hospital Seed Germinator Double Chamber India

Technical Specification :

Each chamber size: 555 (W) x 910 (H) x 605 (D) mm.

For cooling (Low temp.) 5° C to 50° C ± 1° C a sealed compressor unit has been provided to work independently below the chamber.

2A: Heating & Cooling in separate chambers.

2B: Heating & Cooling in both chambers.

Optionals:

Digital temp. indicator-cum-controller with computer data logger

Automatic 0-24 hrs

Timer to regulate cyclic illumination conditions

Chart recorder

Automatic voltage stabilizer - 3 KVA (Must)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003