

**Product Name :**  
Seed Germinator Dual Chamber

**Product Code :**  
NLAB-TESTINGLABM39003



**Description :**

Seed Germinator Dual Chamber

**Technical Specification :**

Double walled. Inside made of stainless & outside of painted M. S.  
Temperature range in one chamber is 5°C to 50°C ± 1°C & Ambient to 50°C ± 1°C in another chamber. Temperature is controlled by Electronic Digital Temperature Controller in both chamber. Each chamber had 14 trays (total 28 trays) and separate air circulation system. Fitted with transparent door for test visibility. Each chamber is fitted with stainless steel water reservoir equipped with immersion heaters to provide uniform humidity of 90% to 95% ± 5% (Not adjustable). To work on 220 V 50 cycles AC

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003