

**Product Name :**  
Self Biasing of Transistor Trainer Digital Meters

**Product Code :**  
NLAB-ELECTRICALAB210064



**Description :**

Self Biasing of Transistor Trainer Digital Meters

**Technical Specification :**

SCOPE OF LEARNING:

Self Biasing of BJT

TECHNICAL SPECIFICATIONS:

Digital Meters:

Voltmeter 20VDC

Ammeter 200uA/2mA DC.

Power Supplies:

DC Supply IC Regulated +12V DC, 250mA.

Operated on Mains power 230V, 50Hz +10%

Components are mounted on the panels are:

Transistor SL100

Required Resistors

**SALIENT FEATURES:**

Front panel built with high class insulated Printed Circuit Board sheet with well printed circuits and symbols.

Fuse for Short Circuit protection

Instruction manual.

Connections are brought out through 4mm Colored Sockets.

Patch Cords 4mm.

The trainer is housed in ABS Plastic cabinet.

Size of the trainer set 12"x8"

## **Naugralabequipments**

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003