

Product Name :
Semi Automatic Double Wall Autoclave

Product Code :
NLAB-LABORATORYEM47013



Description :

Semi Automatic Double Wall Autoclave

Technical Specification :

Semi Automatic Double Wall Autoclave -

Features:-

Reduced loading height

µp based PID controller with dual display for Temp. & Time

Unique pedal free design

Fully SS 304 Construction

No Hinge or shaking lids

Specification:-

Capacity
30 x 500mm Dimension
Power Consumption
Operating / Pressure

021001g Temperature
2301gAC
Digital
0351g Weight

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003