

Product Name :
Semi Fowler Bed Electric

Product Code :
NLABMEDICAL-M2790003

Description :

Semi Fowler Bed Electric India, Semi Fowler Bed Electric manufacturer India, suppliers, Medical Semi Fowler Bed Electric India, Hospital Semi Fowler Bed Electric India

Technical Specification :

Overall Size:

(2160 L x 1015 W x 600 H)

Suitable for mattress size:

(L 1980 x W 910)

Two section top of uniformly perforated CRCA sheet.

Backrest section maneuvered by the help of Actuators to achieve smooth and efficient operation without any jerks

ABS bows easily removable for access to the patient

The bed is fitted with 125 mm die swiveling castors with center locking

Frame Work made of rectangular M.S. tubes

Overall approx size: 1900mm x 900mm x 600mm

Perforated CRC sheet top in two sections

Backrest positions obtained by electric actuator operator by hand remote system.

Telescopic IV Rod with Four Locations

S.S collapsible safety side railings

Bed mounted on 100mm deluxe castors, two with brakes

ABS molded head and foot bows

Pre Treated & Epoxy Powder Coated Finish

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003