

Product Name : Semi Fowler Bed Super	Product Code : NLABMEDICAL-M2800004
Description : Semi Fowler Bed Super India, Semi Fowler Bed Super manufacturer India, suppliers, Medical Semi Fowler Bed Super India, Hospital Semi Fowler Bed Super India	
Technical Specification : Two section top of uniformly perforated CRCA sheet Backrest section maneuvered by the help of Actuators to achieve smooth and efficient operation without any jerks ABS bows easily removable for access to the patient The bed is fitted with 125 mm die swiveling castors with center locking Frame Work made of rectangular M.S. tubes Perforated CRC sheet top in two sections Backrest positions obtained by electric actuator operator by hand remote system Telescopic IV Rod with Four Locations S.S collapsible safety side railings Bed mounted on 100mm deluxe castors, two with brakes ABS molded head and foot bows Pre Treated & Epoxy Powder Coated Finish approximate Size: mattress size {2160(L)× 1015(W) ×600(H)} mm: {1980(L)×910(W)} mm.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003