

Product Name :
Semi Fowler Bed Super

Product Code :
NLABMEDICAL-M2800004

Description :

Semi Fowler Bed Super India, Semi Fowler Bed Super manufacturer India, suppliers, Medical Semi Fowler Bed Super India, Hospital Semi Fowler Bed Super India

Technical Specification :

Two section top of uniformly perforated CRCA sheet

Backrest section maneuvered by the help of Actuators to achieve smooth and efficient operation without any jerks

ABS bows easily removable for access to the patient

The bed is fitted with 125 mm die swiveling castors with center locking

Frame Work made of rectangular M.S. tubes

Perforated CRC sheet top in two sections

Backrest positions obtained by electric actuator operator by hand remote system

Telescopic IV Rod with Four Locations

S.S collapsible safety side railings

Bed mounted on 100mm deluxe castors, two with brakes

ABS molded head and foot bows

Pre Treated & Epoxy Powder Coated Finish

approximate Size: mattress size

{2160(L)× 1015(W) ×600(H)} mm: {1980(L)×910(W)} mm.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003