

Product Name :
Senior Rotary Microtome

Product Code :
NLAB-LABORATORYEM28030



Description :

Senior Rotary Microtome

Technical Specification :

Vibration free operation.

Smooth Drive Wheel with Locking.

Conveniently located feed indicator.

Automatic feed release mechanism.

Sturdy, rigid construction.

Technical Specifications :

Section thickness range: 1-50 microns, 1-50 microns, 1-20 microns

Section thickness setting: In steps 1 micron

Knife length: 120mm

Operation: Manual

Block (Specimen) Holder: Set of three

Object feed: With wheel

Vertical specimen stroke: 55mm, 45mm

Max. specimen size: 32mm, 40x55mm

Adjustment with locks: Provided

Thickness control: Manual

Specimen orientation: 360°

Knife angle: 0 to 10°

Razor holder grip: Single, Double, Double

Dimensions: 300x300x240 mm

Weight 35 kgs: 30 kgs

Dust cover: Provided

Instruction manual: Provided

L Moulds: Optional 50x25x20mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003