

**Product Name :**  
Series Ice Maker

**Product Code :**  
NLAB-SCIENTIFICEQ33002



**Description :**

Series Ice Maker

**Technical Specification :**

series Ice Maker

- \* Use the imported Speed Reducer from Haitec in Italy and Motor from GGM in Korean.
- \* Evaporator are made of Copper.
- \* Compact design , Luxurious appearance.
- \* Full automatic running and easy to use.
- \* A power-overloading protecting device, water shortage auto-detecting, system and ice storage full control device.
- \* High efficiency CFC free compressor, economical and low energy consumption.

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com  
**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

## Naugralabequipments