

Product Name :
Series Ice Maker

Product Code :
NLAB-SCIENTIFICEQ33002



Description :

Series Ice Maker

Technical Specification :

series Ice Maker

- * Use the imported Speed Reducer from Haitec in Italy and Motor from GGM in Korean.
- * Evaporator are made of Copper.
- * Compact design , Luxurious appearance.
- * Full automatic running and easy to use.
- * A power-overloading protecting device, water shortage auto-detecting, system and ice storage full control device.
- * High efficiency CFC free compressor, economical and low energy consumption.

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments