

Product Name :
Series Liquid Manometer

Product Code :
NLAB-TECHNICALAB26004



Description :

Series Liquid Manometer

Technical Specification :

A range of general purpose laboratory manometers using liquid displacement to measure differential pressure.

Capabilities:

Low cost and easy to use instruments

Wide range of pressures accommodated using different manometer fluids

A range of manometers measuring differential water pressures up to approximately 12.5m H₂O. Scales are graduated in 1mm divisions.

1 metre scale differential open water manometer

1 metre scale differential pressurised water manometer (the air space above the tubes can be pressurised with the pump supplied)

1 metre scale differential water over Mercury manometer

500mm scale differential water over Mercury manometer

500mm scale differential kerosene over water manometer

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

