


Product Name : Set of Resonance	Product Code : NLAB-SCHOOLABS15004
	
Description :	
Set of Resonance	
Technical Specification :	
<p>The trial consists of. Resonance tube with scale length 100 cm, water cylinder. Made of good quality aluminum, metal base and pole with adjustable handle. It is used to test the resonance of sound waves in one closed end tube. The length of the test tube can be adjusted by adjusting the water level.</p> <p>Accessories (additional purchase)</p> <p>Tuning Fork</p> <p>Tuning</p> <p>Fork</p>	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

