

Product Name :
Shaking Water Bath with Universal Spring Platform

Product Code :
NLAB-EDUSCHOOLEQ45021



Description :

Shaking Water Bath with Universal Spring Platform

Technical Specification :

This compact shaking water bath utilizes the well proven, belt-less motor drive system from our popular Orbi Shaker Series, with the addition of precision temperature control up to 80°C.

A hinged and removable lid is also included to reduce energy waste, prevent splashing and divert condensation away from samples.

The seamless, stainless steel chamber has a 12 liter capacity and includes a universal spring-loaded platform for securing 15 & 50ml test tubes, as well as Erlenmeyer flasks and media bottles up to 1 liter.

Optional test tube racks are available for microtubes and custom-designed racks can also be made available by request.

Easy to use controls allow the user to digitally select and monitor time, temperature and shaking speed, while a constant monitoring system ensures accuracy over each parameter.

--	--

Ambient Range: 80°C	
Speed Range:	
0.5 bitm / 1 in.	
Temp. Range (air): 37°C	
Temp. Increments:	
Temp. Range (air): 1°C	
Capacity:	
Dimensions (W x H x D): 31 x 34.5 x 30 cm)	
Material: Steel	

Chamber Dimensions	24 x 15.25 cm
Timer	to 48 hours, or continuous
Weight	6.5 kg
Electrical	230V, 50-60Hz, 400W

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003