

**Product Name :**  
Shear Fatigue Tester Machines

**Product Code :**  
NLAB-TESTINGLABM66001



**Description :**

Shear Fatigue Tester Machines

**Technical Specification :**

Shear Fatigue Tester. The provided tester offers a low-cost method to check the peel and shear strength of fasteners. Manufactured using quality-assured components and modern methodology at our in-house production unit,

**Features:**

Portable and compact

High sensitivity

Accurate result

Rugged design

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003