

Product Name :
Shearing Device

Product Code :
NLAB-EDUCATIONTVET13026



Description :

Shearing Device

Technical Specification :

microprocessor controlled, touch screen display, feed speed 0.00001 to 15.0000 mm / min

10: 1 reducer

Comparator 25 mm x 0.01 mm / division

Comparator 10 mm x 0.01 mm / division

Torque ring 3 kN

Set of masses to apply a vertical stress

50 mm sample:

Shear box Ø50 mm

Cutting bag Ø50 mm

Extruder Ø50 mm

Spare porous stone Ø50 mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003