

Product Name : Shop Microscope	Product Code : NLAB-EDUSCHOOLEQ1730128
Description :	
Shop Microscope	
Technical Specification :	
<p>This small hand-held shop microscope is available with 20x, 40x, 50x or 100x magnification and includes a built-in eyepiece reticle for measuring as well as a penlight in a swivel mount for illumination.</p> <p>Weighing less than 1.5 pounds and measuring 9 inches tall, these lightweight shop microscopes are great for mobile use when protected by the accompanying carrying case.</p> <p>Shop microscopes are frequently used in the textile industry and in industrial settings for making quick measurements and quality control inspections.</p> <p>Each microscope includes its own sturdy carrying case.</p>	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003