

**Product Name :**  
Shop Microscope

**Product Code :**  
NLAB-EDUSCHOOLEQ1730128

**Description :**

Shop Microscope

**Technical Specification :**

This small hand-held shop microscope is available with 20x, 40x, 50x or 100x magnification and includes a built-in eyepiece reticle for measuring as well as a penlight in a swivel mount for illumination.

Weighing less than 1.5 pounds and measuring 9 inches tall, these lightweight shop microscopes are great for mobile use when protected by the accompanying carrying case.

Shop microscopes are frequently used in the textile industry and in industrial settings for making quick measurements and quality control inspections.

Each microscope includes its own sturdy carrying case.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003