

Product Name : Siemens 300 Ventilator	Product Code : NLABMEDICAL-M3390002
---	---

--

Description :

Siemens 300 Ventilator India, Siemens 300 Ventilator manufacturer India, suppliers, Medical Siemens 300 Ventilator India, Hospital Siemens 300 Ventilator India

Technical Specification :

Auto mode.
Neonate to adult patient range.
Wide range of ventilation modes.
Easy to use.
O2 air gas module filters bacteria.

Naugralabequipments
Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003