

Product Name :
Siemens 900C Ventilator

Product Code :
NLABMEDICAL-M3390003

Description :

Siemens 900C Ventilator

Technical Specification :

Gas consumption: Approximately equal to preset aspiratory minute volume or equal to flow meter setting.

Gas supply via high pressure inlet: From central gas supply, compressor or gas cylinder.

Electronic gas supply unit, inlet pressure: 2.5-5 bar (250-500 kPa, 35-70 psi).

Other gas supply inlet pressure: 2.5-7 bar (250-700 kpa, 35-100 psi).

Gas supply via low pressure inlet: Vial flow meters.

Power supply: (100), 110, 120, 220, 240, 50 or 60 hz.

Power consumption: Approximately 40 W.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003