

Product Name :
Sieve Shaker

Product Code :
NLAB-TESTINGLABM109071



Description :

Sieve Shaker

Technical Specification :

Sieve Shaker

Sieve Shakers, hand operated or motorized, have been designed to provide a standard sieving action for perfect and quick sieving of soil and other engineering materials. The motorized unit save considerable time. Time switch adjustable from 0-60 min., in time intervals of 5 min is provided as standard along with the motorized sieve shakers.

Sieve Shaker (Manually) 20cm

This is a light, portable but sturdy sieve shaker suitable for bench mounting the side to side movement to the carrier which can take upto 7 sieves of 200mm dia is through a train of gears operated by a hand wheel. A heavy fly wheels ensures smooth operation.

Sieve Shaker (Manually) 20cm, 30cm

Same as above but carrier table & holder suitable to take 20 cm & 30cm (200mm & 300mm) dia sieves.

Sieve Shaker (Manually) 20cm & 30cm

Sieves Box Type

Same as above but carrier table & holder suitable to take 20 cm & 30cm (200mm & 300mm) dia sieves.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003