

**Product Name :**  
Sieve Shaker

**Product Code :**  
NLAB-TESTINGLABM109071



**Description :**

**Sieve Shaker**

**Technical Specification :**

**Sieve Shaker**

Sieve Shakers, hand operated or motorized, have been designed to provide a standard sieving action for perfect and quick sieving of soil and other engineering materials. The motorized unit save considerable time. Time switch adjustable from 0-60 min., in time intervals of 5 min is provided as standard along with the motorized sieve shakers.

**Sieve Shaker (Manually) 20cm**

This is a light, portable but sturdy sieve shaker suitable for bench mounting the side to side movement to the carrier which can take upto 7 sieves of 200mm dia is through a train of gears operated by a hand wheel. A heavy fly wheels ensures smooth operation.

**Sieve Shaker (Manually) 20cm, 30cm**

Same as above but carrier table & holder suitable to take 20 cm & 30cm (200mm & 300mm) dia sieves.

**Sieve Shaker (Manually) 20cm & 30cm**

**Sieves Box Type**

Same as above but carrier table & holder suitable to take 20 cm & 30cm (200mm & 300mm) dia sieves.

## **Naugralabequipments**

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003