

**Product Name :**  
Sieve With Three-Dimensional Movement, Type AS 200 Basic

**Product Code :**  
NLAB-EDUCATIONTVET13014



**Description :**

Sieve With Three-Dimensional Movement, Type AS 200 Basic

**Technical Specification :**

Measuring range: 20 µm to 25 mm (depending on sieve and type of sample)

Maximum load: 3 kg

Capacity: 8 H50 mm sieves or 16 H25 mm sieves

Setting the sieving parameters:

Amplitude: digital, 1 to 100% (~ 3 mm)

Timer: digital, 1 to 99 minutes

Usable screen diameters: 100 to 203 mm

Height of the screening column: max. 450 mm

Dimensions (WxHxD): 417 x 212 x 384 mm

Weight: ~ 35 kg

Power supply: 230V - 50Hz

"Economy" fastening system, for dry sieving, for 100/150/200/203 mm sieves

Collection bottom, diameter 200 mm, height 50 mm

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003