

Product Name :
Silicon Tubing Rubber

Product Code :
NLAB-SCIENTIFICEQ31004



Description :

Silicon Tubing Rubber

Technical Specification :

Highly elastic, non-toxic, heat resistant from -50°C to 250°C , repeatedly autoclave, flexible odorless and transparency of the tubes allows clear visibility of the flowing fluid.

Sizes (bore $\times$ wall thickness) available are:-

2 $\times$ 1.0

2 $\times$ 2.0

3 $\times$ 1.0

4 $\times$ 1.5

5 $\times$ 1.5

6 $\times$ 1.5

7 $\times$ 2.0

8 $\times$ 1.5

8×2.0

10×1.5

10×2.02

14×3.0

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003