

Product Name :
Simple Galvanoscope

Product Code :
NLAB-EDUSCHOOLEQ1740119



Description :

Simple Galvanoscope

Technical Specification :

Simple Galvanoscope for Physics Lab comprises wooden former with terminals & engraved metal circle marked 0 – 180 – 0, divided & figured every 45°, with magnetic pivot, wound with 30 turns enamelled copper wire & complete with 75 mm magnetic needle with brass bearing. Maximum opening voltage 4.5 V.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003