

Product Name :
Simple Lever

Product Code :
NLAB-LABORATORYEM77023

Description :

Simple Lever

Technical Specification :

Stout construction duly chromium plated frame with stop at one end, interchangeable " stem.

The lever fits through axle pivoted between centre.

A thumb screw on the axle permits the lever to be adjusted for length.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003