

Product Name : Simple Lever	Product Code : NLAB-LABORATORYEM77023
Description :	
Simple Lever	
Technical Specification :	
Stout construction duly chromium plated frame with stop at one end, interchangeable " stem.	
The lever fits through axle pivoted between centre.	
A thumb screw on the axle permits the lever to be adjusted for length.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003