

**Product Name :**  
Single Cylinder Two Stage Air Compressor Test Rig

**Product Code :**  
NLAB-TECHNICALAB25008



**Description :**

Single Cylinder Two Stage Air Compressor Test Rig

**Technical Specification :**

One experimental type 2 stage reciprocating air compressor, air cooler, suitable for a pressure of 12 Kg/Cm<sup>2</sup> having a displacement of 0.4 cubimetre per minute, driven by squirrel cage induction motor 5 HP turnion mounted electric dynamometer type V belt drive and fitted with torque arm and spring balance.

The above motor suited 400/440v 3-phase 50 Cycles A.C. Supply.

The motor is provided with switch foundation bolted.

Accessories:

The Equipment complete with the following standard accessories comprising the following: Two Pressure gauges for measuring the pressure of each stage.

One shut off valves, pressure released valve.

A receiver fitted with oil and water separator a drain cock and automatic pressure switch.

Air Intake Measuring Arrangement.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

