

Product Name :
Site Dumper

Product Code :
NLAB-VOCATIONALTVET15005



Description :

Site Dumper

Technical Specification :

Payload capacity 3000kg.

Skip capacity (litre) 1500.

Chassis: Articulated Frame.

Gearbox: Heavy duty, 5+5 gearbox.

Bucket tipping system: Hydraulic tipping hopper.

Diesel engine, 40hp, water cooling, 30kw/2300rpm.

Driving 4x4 WD.

Tires Four wheels 12.5-16.

Turning radius (m) 4.5.

Hydraulic steering.

Weight (kg) 2250.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003