

Product Name :
Skid Resistance and Friction Tester

Product Code :
NLAB-AGRETLAB160021



Description :

Skid Resistance and Friction Tester

Technical Specification :

Used for the measurement of surface friction properties, the apparatus is suitable for both site and laboratory applications and for polished stone value tests using curved specimens from accelerated polishing tests. The equipment is supplied with:

Skid Resistance And Friction Tester

Additional scale for tests on polished stone value specimens

6 rubber sliders for site use, complete with conformity certificate

Thermometer -10 to +110°C for surface temperature measurement

1 litre washing bottle, for surface wetting

Tool set with case, for machine assembly

Rule for sliding length verification

Rocker weight : 1500 +/- 30 g

Distance between rocking centre & center of gravity : 410 +/- 5 mm

Positive static pressure on pavement : 22.2 +/-0.5 kn

6 rubber slider and base plate dimension (mm) : 700x360x70

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003