

Product Name :
Slit Lamp

Product Code :
NLABMEDICAL-M3460002

Description :

Slit Lamp

Technical Specification :

Magnification: 5x,10x,15x,20x & 25x continuous zoom

Interpupillary Distance: 52mm to 80mm

Camera: CCD Digital Camera

Slit Width Adjustment: 0 to 12mm without steps

Slit Length Adjustment: 0 to 12mm continuous

Diameter of Diaphragm: 12mm, 10mm, 8mm, 5mm, 3mm, 2mm, 1mm and 0.2mm

Angle of Rotation of Slit: $\pm 90^\circ$

Tilt of Slit: Horizontal in steps of $0^\circ, 5^\circ, 10^\circ, 15^\circ$, Vertical in steps $0^\circ, 5^\circ, 10^\circ, 15^\circ, 20^\circ$

Filters: Cobalt Blue, Green, Neutral Density, Heat absorbing

Light Source: Halogen Lamp 12v-50w

Intensity Control of Illumination: Low, Medium and High

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003