

Product Name :
Slotted Brass Weight Set

Product Code :
NLAB-EDUSCHOOLEQ1740120



Description :

Slotted Brass Weight Set

Technical Specification :

Slotted Brass Weight Set – Made of brass. This set includes nine slotted masses each of 10 g with a 10 g wire hanger making total mass of 100 g; masses do not fall off when assembled in a set. Tolerance $\pm 3\%$. Also available in 200g, 500g, 1000g.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003