

**Product Name :**  
Smoke Machine Evaporative System Tester

**Product Code :**  
TK554



**Description :**

Smoke Machine Evaporative System Tester

**Technical Specification :**

Approved for testing vehicle`s evaporative emissions (EVAP) System.  
Two-minute Pass / Fail feature verifies .040" or .020" leaks in OBD-II vehicles.  
Smoke out OBD-II EVAP codes.  
All accessories included.  
Patented Ultra Trace UVT technology deposits a fluorescent dye at the exact location of the leaks.  
500 tests per refill.  
User refillable.  
A tool for OBD-II EVAP testing & general-purpose leak testing.  
Nitrogen-based for safe EVAP testing & uses shop-air for general purpose leak testing.  
Microprocessor technology.  
Dimensions: 15" x 8" x 15"  
Shipping Weight: 25 lb.  
Power Supply: 12 volts DC  
Power Consumption: 15 amps  
Supply Pressure: 13 in./H2O  
Supply Volume: 15 LPM  
Operating Temp. Range: 45Â° to 140Â°  
Smoke Supply Line: 10`  
Power Supply Line: 10`  
Remote Starter Cable: 10`

The Vacuity is heavy-duty all metal construction and has been designed to deliver years of trouble-free operation with a yearly smoke solution change as its only required maintenance.

With full accessories:

All components and leads made to protect the safety of the students according to the laws in force.

It provided with manual in English language

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003