


Product Name :
Soil Moisture Meter Testing

Product Code :
NLAB-TESTINGLABM30001



Description :

Soil Moisture Meter Testing

Technical Specification :

This is a simple water sensor, can be used to detect soil moisture

Module Output is high level when the soil moisture deficit, or output is low

Can be used in module plant waterer device, and the plants in your garden no need people to manage

Operating voltage: 3.3V~5V

Dual output mode, analog output more accurate

A fixed bolt hole for easy installation

With power indicator (red) and digital switching output indicator (green)

Having comparator chip, stable

Panel PCB Dimension: Approx. 3cm x 1.5cm

Soil Probe Dimension: Approx. 6cm x 3cm

Cable Length: Approx.21cm

Interface Description(4-wire):

VCC: 3.3V-5V

GND: GND

DO: digital output interface(0 and 1)

AO: analog output interface

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003