

Product Name : Soil Moisture Suction Sand Table	Product Code : NLAB-TECHNICALAB27005
Description :	
Soil Moisture Suction Sand Table	
Technical Specification :	
<p>Sand table, 300mm diameter Sliding suction control chamber giving tension range 0-800mm of water (using sandy-loam) Complete with three soil sample retaining rings (80mm diameter), instructions and manual. These equipment consists of a free-standing circular tank filled with saturated sand (not supplied). Equilibrated soil samples placed on the sand are then subjected to a known variable suction range of 0- 800mm of water. The suction is applied by means of a buried porous pipe and an associated sliding control chamber. In order to maintain the suction an uninterrupted low volume water supply is required. A lightweight cover is provided to prevent evaporative water loss during experiments and the table is large enough to permit the simultaneous evaluation of several samples.</p>	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003