

Product Name :
Soil Moisture Tester**Product Code :**
SOLE0003**Description :**

Soil Moisture Tester

Technical Specification :

Soil Moisture Tester

Features :

The speedy Soil Moisture Tester measures specimen weight before and after drying, then get percentage of water in the specimen

Technical Specifications :

Capacity : 110 g

Division : 5 mg

Repeatability : 0.2%

Readability : 0.02%

Temperature Range : Room Temp. to 200°C

Temperature accuracy : $\pm 1^\circ\text{C}$

Lamp : Halogen

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-0171-2643080, 2601773