

**Product Name :**  
Soil Testing Advanced Dynamic Triaxial Testing System

**Product Code :**  
NLAB-EDUSCHOOLEQ200001



**Description :**

Soil Testing Advanced Dynamic Triaxial Testing System

**Technical Specification :**

Axial force and axial deformation are applied through the base of the cell.

The Advanced Dynamic Triaxial Testing System is a high-end, no compromise testing apparatus combining a triaxial cell with a dynamic actuator capable of applying load, deformation and stresses up to 10Hz.

The cell itself is screw-driven from an integral base unit housing the motor drive.

The system can be combined with a dynamic cell pressure actuator such that cell pressures may also be applied dynamically up to fundamental frequency of the machine (i.e. 2Hz, 5Hz or 10Hz).

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003