

Product Name :
Soil Testing Extractor Frame Hydraulic

Product Code :
NLAB-EDUSCHOOLEQ200011



Description :

Soil Testing Extractor Frame Hydraulic

Technical Specification :

We have in our range Extractor Frame – Hydraulic that is highly acknowledged in the market for their premium quality and superior functionality. These frames are comprises 5 ton capacity hand operated hydraulic jack that is mounted on a suitable frame two plungers. One is of 100 mm dia and the other of 150mm dia moulds are provided on one thrust plate for 150 mm dia. The height of thrust plate is adjustable and a set of plunger adaptor thrust plates for 38mm, 50mm and 75mm is incorporated.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003