

Product Name :
Soil Testing Load Frame

Product Code :
NLAB-EDUSCHOOLEQ200020



Description :

Soil Testing Load Frame

Technical Specification :

Motorised cum hand: Tow pillar type

4 speeds: Ideal for use in mobile laboratory as it is light weight & sturdy. This Load frame is designed for conducting Unconfined, Triaxial CBR and other routine tests.

Loading system comprises of a screw jack with detachable handle. The lower platen moves up and down. The top bracket is adjustable for vertical clearance and has an adapter for connecting standard Proving ring.

A dial gauge mounting bracket is provided on one to the two pillars. Rate of strain 1.25, 1.50 & 2.50 mm/min. Suitable for operation on 220V, 1 Phase 50 cycles, AC supply.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003