

Product Name :
Soil Testing Marsh Cone Funnel

Product Code :
NLAB-EDUSCHOOLEQ200021



Description :

Soil Testing Marsh Cone Funnel

Technical Specification :

This cone is used to find out viscosity of bentonite slurry and like material. The marsh cone is 6 inch in diameter at the top and 11 inch(280mm) long, and tapers to join a tube 2 inch long and 3/16 inch side dia. The capacity of the funnel is 1500 cc. Time in seconds required to flow out 1000 cc of slurry from cone is measured as funnel viscosity of the material.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003