

Product Name :
Soiled Linen Trolley Canvas Bag Round

Product Code :
NLABMEDICAL-M2620013



Description :

Soiled Linen Trolley Canvas Bag Round

Technical Specification :

Specification :-

Frame is made of Round tube.

It comes with good quality superior canvas bag.

The complete soiled linen trolley is mounted on 5 cms castors.

We are Manufacturer, Exporters and Suppliers of Soiled Linen Trolley Canvas Bag Round at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003