

Product Name :
Soiled Linen Trolley Plastic Bucket

Product Code :
NLABMEDICAL-M2620015



Description :

Soiled Linen Trolley Plastic Bucket

Technical Specification :

Specification :-

Frame is made of round tube.

It comes with good quality superior plastic bucket.

Push handle is available at one end.

The complete soiled linen trolley is mounted on 10 cms castors.

We are Manufacturer, Exporters and Suppliers of Soiled Linen Trolley Plastic Bucket at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

