

Product Name :
Sole Adhesion Tester Machines

Product Code :
NLAB-TESTINGLABM96001



Description :

Sole Adhesion Tester Machines

Technical Specification :

Sole Adhesion Tester technology support and finds use in meeting the measuring adhesion strength of soles at toe and heel area. Finding extensive use in shoe manufacturing units, these testers provide for proven performance and are also used in laboratory conditions. Some of its features include maximum force reading of dial indicator is and very reliable tester functioning, supports demands of testing strength of adhesion of sole to its upper; accurate measuring standards

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003