

Product Name :
Solubility Product With Cobra 4

Product Code :
NLAB-EDUCATIONTVET12003

Description :

Solubility Product With Cobra 4

Technical Specification :

Each experience includes:

Immersion thermostat

Cobra4 mobile-link 2 - mobile recorder with SD card, rechargeable batteries, USB cable and software

Cobra4 conductivity / temperature sensor (Pt1000)

Conductivity-temperature probe

Magnetic agitator

Thermostat tank

Bunsen stand, 75 cm

Universal pliers

Mortar with pestle

Double nuts

Narrow neck Erlenmeyer flask 100 ml

Spoon with spatula

Magnetic rod 30 mm

500 ml wash bottle

Powder funnel

Recent configuration laptop

Consumables:

Standard conductivity solution, 1.413 $\mu\text{S} / \text{cm}$, 25 sachets of 20ml

Calcium carbonate, 500 g

Calcium fluoride, 100 g

Weighing platform, 84x84mm, 25 Pcs

Balance 400g / 1mg (Cfr item H-2)

Distiller (Cfr item H-4)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003