

**Product Name :**  
Solvent Recovery Unit

**Product Code :**  
NLAB-BTMENLAB150022



**Description :**

Solvent Recovery Unit

**Technical Specification :**

Used to recover the solvent liquid after its use for the extraction tests. This unit has been designed to recover non-flammable solvents and consists of two stainless steel chambers, one for dirty solvent and the other for the cleaned solvent. An electric heater in the left-hand chamber distils the solvent, which then passes through a water cooling system and drops into the second chamber ready for re-use in a new test. Once the process is completed, a temperature switch automatically stops the heating elements. Supplied complete with 10 m plastic tube, tube clamps, sieve insert 0.6 mm opening and one lid.

Max. temperature: 150°C

Power: 1200 W

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

