

Product Name :
Sonometer Mini

Product Code :
NLAB-SCIENTIFICEQ26002



Description :

Sonometer Mini

Technical Specification :

Compact in size with overall 50c.m scale to maintain a gap between two bridges.

A pulley is fitted at the extreme end of the sonometer to hang the weight set on the sample wire.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003