

Product Name : Soya Peptone	Product Code : NLAB-CHEMICALSIND20019
Description :	
Soya Peptone	
Technical Specification :	
Soya Peptone Type I	
Soya peptone type I, is a soluble end product of the enzymic digestion of soya Peptone meal by papain. It is use as growth stimulant in culture media for the mass cultivation of variety of microorganisms including bacteria and fungi. It is excellent source of mixed peptides, free amino acids and growth factors.	
Quality Control	
pH: 6.0 - 7.0	
Reaction: Reaction of 2% w/v aqueous solution at 25°C.	
Appearance: Yellow to yellowish brown homogenous free flowing characteristic odour of protein, derived from vegetable source.	
Solubility: Freely soluble in distilled/purified water, insoluble in chloroform.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003