

Product Name :
Speciman Jar

Product Code :
NLAB-EDUSCHOOLEQ129001



Description :

Speciman Jar

Technical Specification :

These transparent specimen jars, used for preserving different specimens in the lab are provided with a cap that fits in tight.

These jars provided excellent protection to the species preserved in it.

Speciman Jar ? 50 mm X h 150 mm.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003